

Mounting of the mattress


RotoBed®Care comes with a mattress consisting of a main part and a foot part.

The main mattress has 5 buttons on the back. First place the three lower buttons in the grooves at the bottom of the bed. Slide the mattress towards the foot end. Now press the two buttons in the centre of the mattress into the corresponding holes in the middle of the bed.

Similarly, the foot mattress has two buttons, which fit into the two round holes in the bed.

Mounting of side rails

RotoBed® care beds can be delivered with side rails made of metal or oak laminate. Regardless of the material the side rails are mounted and removed in the same way:



Press the small fitting underneath while lifting the bed rail.

The side rails can simply be put back into place and will automatically lock.

Castors and brakes

As standard, the bed comes with integral Tente Linea double castors. The castors are rotatable. The brakes work in pairs and can be operated from both sides of the bed.



Note: RotoBed®Care should not be used by people without an introduction in the use of the bed.

The bed is tested according to the latest applicable standards and rules for care beds. Be sure that neither children nor animals are under the bed when the bed's parts are set. When the bed's leg part is raised, do not overload it with heavy objects or persons.

The bed must be used by users within the following criteria: Weight 40-190 kg. Height \geq 146 cm. BMI \geq 17.

Quick guide for RotoBed® Care

Exit: Push the **Chair Position** key until the bed stops in comfort position. Release the rotation mechanism by slightly pushing the grip with your foot. Now slowly turn the bed until the grip automatically locks. Push **Chair Position** key again until the exit position. To facilitate the exit you can, by means of the remote control, move the backrest the last bit forward and raise the bed a little in height.




Entry: Help the user getting seated as close to the back as possible. Push the **Sleeping Position** key until the bed stops in comfort position. Release the rotation mechanism by slightly pushing the grip by foot. Now slowly turn the bed until the grip automatically locks. Push the **Sleeping Position** key again until the bed is rotated and horizontal.



Standard seat height without mattress is 32 cm above the floor. For individual settings, see user manual.

Remote control



RotoBed®Care can turn either right or left. Therefore, the remote control can also be connected to both the right and left sides of the bed. Optional remote control can be purchased as desired. When the remote control is not in use, place it on the side rail.